

## Single tooth restorations with Astra Tech Implant System™

The first studies on Astra Tech Implant System™ evaluating the outcome of single tooth rehabilitations were published already in 1997<sup>1-4</sup>. High implant survival rates (more than 97.8%), extremely small changes in the marginal bone levels (range: 0.0–0.3 mm), and surprisingly small numbers of abutment screw loosening were reported after 1–2 years in function.

Today, there are over 40 articles, including 10 year results<sup>5</sup>, published in per-reviewed journals reporting on the proof of concept using Astra Tech implants for single tooth rehabilitations. Examples are:

- clinical performance in posterior<sup>1,6-8</sup> and or anterior<sup>2,4,9-19</sup> locations
- clinical performance when applying 1-stage surgery<sup>6,8-10,12,20</sup> combined with immediate loading<sup>11,15,21-25</sup> or after placement in extraction sockets<sup>12,20,24,26-30</sup>
- meticulous esthetical evaluations<sup>20,23,25,31-34</sup>
- the long-term outcome<sup>5,8,12,14,16,35-37</sup>
- bone regeneration using different sinus lift techniques<sup>6,7,38-40</sup>
- periodontally compromised<sup>41-43</sup> or diabetic patients<sup>44</sup>
- microbiological<sup>45,46</sup> or antibacterial<sup>47</sup> focuses
- patient evaluated satisfaction or quality of life<sup>16,31,48,49</sup>

Implant survival rates are in general 94–100%, irrespective of surgical technique or location in the mouth. The documentation clearly shows that the supporting bone is maintained around Astra Tech single tooth implants. Maxillary anterior implants can actually obtain a gain in marginal bone levels over a 5 years period<sup>16</sup>. Hereof, esthetics is also well maintained and rated high by patients and dentists.

1. Karlsson U, Gotfredsen K, Olsson C. Single-tooth replacement by osseointegrated Astra Tech dental implants: a 2-year report. *Int J Prosthodont* 1997;10(4):318-24. (ID No. 75067) [Abstract in PubMed](#)
2. Kempainen P, Eskola S, Ylipaavalniemi P. A comparative prospective clinical study of two single-tooth implants: a preliminary report of 102 implants. *J Prosthet Dent* 1997;77(4):382-7. [Abstract in PubMed](#)
3. Norton MR. The Astra Tech single-tooth implant system: a report on 27 consecutively placed and restored implants. *Int J Periodontics Rest Dent* 1997;17(6):575-83.
4. Palmer RM, Smith BJ, Palmer PJ, Floyd PD. A prospective study of Astra single tooth implants. *Clin Oral Implants Res* 1997;8(3):173-9. (ID No. 75182) [Abstract in PubMed](#)
5. Gotfredsen K. A 10-year prospective study of single tooth implants placed in the anterior maxilla. *Clin Impl Dent Rel Res* 2009;E-pub Aug 6, DOI: 10.1111/j.1708-8208.2009.00231.x. [Abstract in PubMed](#)
6. Diss A, Dohan DM, Mouhyi J, Mahler P. Osteotome sinus floor elevation using Choukroun's platelet-rich fibrin as grafting material: a 1-year prospective pilot study with microthreaded implants. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2008;105(5):572-9. [Abstract in PubMed](#)
7. Fermergård R, Åstrand P. Osteotome sinus floor elevation and simultaneous placement of implants—a 1-year retrospective study with Astra Tech implants. *Clin Impl Dent Rel Res* 2008;10(1):62-9. [Abstract in PubMed](#)
8. Norton MR. Multiple single-tooth implant restorations in the posterior jaws: maintenance of marginal bone levels with reference to the implant-abutment microgap. *Int J Oral Maxillofac Implants* 2006;21(5):777-84. (ID No. 78773) [Abstract in PubMed](#)
9. Cooper L, Felton DA, Kugelberg CF, Ellner S, Chaffee N, Molina AL, et al. A multicenter 12-month evaluation of single-tooth implants restored 3 weeks after 1-stage surgery. *Int J Oral Maxillofac Implants* 2001;16(2):182-92. (ID No. 75410) [Abstract in PubMed](#)
10. Cooper LF, Ellner S, Moriarty J, Felton DA, Paquette D, Molina A, et al. Three-year evaluation of single-tooth implants restored 3 weeks after 1-stage surgery. *Int J Oral Maxillofac Implants* 2007;22(5):791-800. (ID No. 78988) [Abstract in PubMed](#)
11. De Kok JJ, Chang SS, Moriarty JD, Cooper LE. A retrospective analysis of peri-implant tissue responses at immediate load/provisionalized microthreaded implants. *Int J Oral Maxillofac Implants* 2006;21(3):405-12. (ID No. 78727) [Abstract in PubMed](#)
12. Gotfredsen K. A 5-year prospective study of single-tooth replacements supported by the Astra Tech implant: a pilot study. *Clin Impl Dent Rel Res* 2004;6(1):1-8. (ID No. 78273) [Abstract in PubMed](#)
13. Norton MR. Marginal bone levels at single tooth implants with a conical fixture design. The influence of surface macro- and microstructure. *Clin Oral Implants Res* 1998;9(2):91-9. [Abstract in PubMed](#)
14. Norton MR. Biologic and mechanical stability of single-tooth implants: 4- to 7-year follow-up. *Clin Impl Dent Rel Res* 2001;3(4):214-20. [Abstract in PubMed](#)
15. Norton MR. A short-term clinical evaluation of immediately restored maxillary TiOblast single-tooth implants. *Int J Oral Maxillofac Implants* 2004;19(2):274-81. (ID No. 78173) [Abstract in PubMed](#)
16. Palmer RM, Palmer PJ, Smith BJ. A 5-year prospective study of Astra single tooth implants. *Clin Oral Implants Res* 2000;11(2):179-82. (ID No. 75352) [Abstract in PubMed](#)
17. Nisapakulorn K, Suphanantachai S, Silkoessak O, Rattanamongkolgul S. Factors affecting soft tissue level around anterior maxillary single-tooth implants. *Clin Oral Implants Res* 2010;E-pub DOI 10.1111/j.1600-0501.2009.01887.x.
18. Martin R. Astra Tech OsseoSpeed™ 3.0S implant. *Inside Dentistry* 2010;6(4):2-4.
19. Norton MR. The Astra Tech Single-Tooth Implant System: a report on 27 consecutively placed and restored implants. *Int J Periodontics Rest Dent* 1997;17(6):574-83. [Abstract in PubMed](#)
20. Lops D, Chiapasco M, Rossi A, Bressan E, Romeo E. Incidence of inter-proximal papilla between a tooth and an adjacent immediate implant placed into a fresh extraction socket: 1-year prospective study. *Clin Oral Implants Res* 2008;19(11):1135-40. (ID No. 79132) [Abstract in PubMed](#)
21. Bilhan H, Sonmez E, Mumcu E, Bilgin T. Immediate loading: three cases with up to 38 months of clinical follow-up. *J Oral Implantol* 2009;35(2):75-81. [Abstract in PubMed](#)
22. Donati M, La Scala V, Billi M, Di Dino B, Torrisi P, Berglund T. Immediate functional loading of implants in single tooth replacement: a prospective clinical multicenter study. *Clin Oral Implants Res* 2008;19(8):740-48. (ID No. 79065) [Abstract in PubMed](#)
23. Harvey BV. Optimizing the esthetic potential of implant restorations through the use of immediate implants with immediate provisionals. *J Periodontol* 2007;78(4):770-6. [Abstract in PubMed](#)
24. Cooper L, Raes F, Reside G, Soliva Garriga J, Ginder Tarrida L, Wilfang J, et al. Comparison of radiographic and clinical outcomes following immediate provisionalization of single tooth dental implants placed in healed alveolar ridges and extraction sockets. *Int J Oral Maxillofac Implants* 2010;25(6):1222-32.
25. Raes F, Cosyn J, Crommelinck E, Coessens P, De Bruyn H. Immediate and conventional single implant treatment in the anterior maxilla: 1-year result of a case series on hard and soft tissue response and aesthetics. *J Clin Periodontol* 2010;accepted for publication Nov 25, 2010, DOI: 10.1111/j.1600-051x.2010.01687.x.
26. Norton MR, Wilson J. Dental implants placed in extraction sites implanted with bioactive glass: human histology and clinical outcome. *Int J Oral Maxillofac Implants* 2002;17(2):249-57. (ID No. 75419) [Abstract in PubMed](#)
27. Misch C. The use of recombinant human bone morphogenetic protein-2 for the repair of extraction socket defects: a technical modification and case series report. *Int J Oral Maxillofac Implants* 2010;25:1246-52.
28. Acocella A, Bertolai R, Sacco R. Modified insertion technique for immediate implant placement into fresh extraction socket in the first maxillary molar sites: a 3-year prospective study. *Implant Dent* 2010;19(3):220-8. [Abstract in PubMed](#)
29. Kahnberg KE. Immediate implant placement in fresh extraction sockets: a clinical report. *Int J Oral Maxillofac Implants* 2009;24(2):282-8. [Abstract in PubMed](#)
30. Valentini P, Abensur D, Albertini JF, Rocchessani M. Immediate provisionalization of single extraction-site implants in the esthetic zone: a clinical evaluation. *Int J Periodontics Rest Dent* 2010;30(1):41-51. [Abstract in PubMed](#)
31. Palmer RM, Farkondeh N, Palmer PJ, Wilson RF. Astra Tech single-tooth implants: an audit of patient satisfaction and soft tissue form. *J Clin Periodontol* 2007;34(7):633-8. (ID No. 78941) [Abstract in PubMed](#)
32. Cody RP. Esthetics in implant dentistry – a case report using the Astra Tech Zir Abutment in a maxillary anterior single tooth reconstruction. *Us Dentistry* 2006:27-28.
33. Lee DW, Huh JK, Park KH, Chai JK, Kim CK, Moon IS. Comparison of interproximal soft tissue height for single implants and contra-lateral natural teeth. *Clin Oral Implants Res* 2009;20(12):1320-25. [Abstract in PubMed](#)
34. Bressan E, Paniz G, Lops D, Corazza B, Romeo E, Favero G. Influence of abutment material on the gingival color of implant-supported all-ceramic restorations: a prospective multicenter study. *Clin Oral Implants Res* 2010;E-pub Nov 11, 2010, DOI: 10.1111/j.1600-0501.2010.02008.x. [Abstract in PubMed](#)
35. Steveling H, Roos J, Rasmusson L. Maxillary implants loaded at 3 months after insertion: results with Astra Tech implants after up to 5 years. *Clin Impl Dent Rel Res* 2001;3(3):120-4. (ID No. 75414) [Abstract in PubMed](#)
36. Wennström JL, Ekestubbe A, Gröndahl K, Karlsson S, Lindhe J. Implant-supported single-tooth restorations: a 5-year prospective study. *J Clin Periodontol* 2005;32(6):567-74. (ID No. 78476) [Abstract in PubMed](#)
37. Chang M, Wennström JL. Longitudinal changes in tooth/single-implant relationship and bone topography: an 8-year retrospective analysis. *Clin Impl Dent Rel Res* 2010;E-pub Feb 11, 2010, DOI: 10.1111/j.1708-8208.2010.00272.x. [Abstract in PubMed](#)
38. Fermergård R, Åstrand P. Osteotome sinus floor elevation without bone grafts - a 3-year retrospective study with Astra Tech implants. *Clin Impl Dent Rel Res* 2009;E-pub Nov 10, DOI: 10.1111/j.1708-8208.2009.00254.x. [Abstract in PubMed](#)
39. Thor A, Senerby L, Hirsch JM, Rasmusson L. Bone formation at the maxillary sinus floor following simultaneous elevation of the mucosal lining and implant installation without graft material: an evaluation of 20 patients treated with 44 Astra Tech implants. *J Oral Maxillofac Surg* 2007;65(Suppl 1):64-72. (ID No. 78929) [Abstract in PubMed](#)
40. Kahnberg KE, Wallstrom M, Rasmusson L. Local sinus lift for single-tooth implant. I. Clinical and radiographic follow-up. *Clin Impl Dent Rel Res* 2009;E-pub Sep 9, DOI: 10.1111/j.1708-8208.2009.00201.x. [Abstract in PubMed](#)
41. Baelum V, Ellegaard B. Implant survival in periodontally compromised patients. *J Periodontol* 2004;75(10):1404-12. [Abstract in PubMed](#)
42. Ellegaard B, Baelum V, Karring T. Implant therapy in periodontally compromised patients. *Clin Oral Implants Res* 1997;8(3):180-8. (ID No. 75060) [Abstract in PubMed](#)
43. Ellegaard B, Kolsen-Petersen J, Baelum V. Implant therapy involving maxillary sinus lift in periodontally compromised patients. *Clin Oral Implants Res* 1997;8(4):305-15. [Abstract in PubMed](#)
44. Turkyilmaz I. One-year clinical outcome of dental implants placed in patients with type 2 diabetes mellitus: a case series. *Implant Dent* 2010;19(4):323-9. [Abstract in PubMed](#)
45. Puchades-Roman L, Palmer RM, Palmer PJ, Howe LC, Ide M, Wilson RF. A clinical, radiographic, and microbiologic comparison of Astra Tech and Brånemark single tooth implants. *Clin Impl Dent Rel Res* 2000;2(2):78-84. (ID No. 75354) [Abstract in PubMed](#)
46. DeAngelo SJ, Kumar PS, Beck FM, Tatakis DN, Leblebicioglu B. Early soft tissue healing around one-stage dental implants: clinical and microbiologic parameters. *J Periodontol* 2007;78(10):1878-86. [Abstract in PubMed](#)
47. Khoury SB, Thomas L, Walters JD, Sheridan JF, Leblebicioglu B. Early Wound Healing Following One-Stage Dental Implant Placement With and Without Antibiotic Prophylaxis: A Pilot Study. *J Periodontol* 2008;79(10):1904-12. [Abstract in PubMed](#)
48. Goshima K, Lexner MO, Thomsen CE, Miura H, Gotfredsen K, Bakke M. Functional aspects of treatment with implant-supported single crowns: a quality control study in subjects with tooth agenesis. *Clin Oral Implants Res* 2010;21(1):108-14. [Abstract in PubMed](#)
49. Van Lierde KM, Corthals P, Browaeyns H, Mussche P, Van Kerckhove E, De Bruyn H. Impact of anterior single-tooth implants on quality of life, articulation and oromyofunctional behaviour: a pilot study. *J Oral Rehabilitation*, DOI: 10.1111/j.1365-2842.2010.02161.x 2010.

